




Stablelab detects infections with 30 times greater sensitivity than a thermometer.¹



Stablelab is a hand held device that measures Serum Amyloid A, a biomarker of infection, and provides results in 10 minutes stall side.



Detect

Detect and measure infection



Monitor

Monitor response to treatment



Screen

Implement protocols to catch problems early

Incorporate Stablelab into your routine clinical exams in cases such as:

Ascending placentitis
ADR
Cellulitis
Colic
Critical care
Diarrhea
EIPH

General infections
Infectious disease
Joint sepsis/flare
New foal exam
Fever of unknown origin
Maladjusted foal
Metritis

NICU
Peritonitis
Pneumonia
Poor performance
Postoperative follow up
Postpartum
Premature foal

Preoperative screening
Respiratory disease
Rhodococcus equi
Rotavirus
Shipping fever
Strangles
Viral infections
Umbilical infection

hello@stablelab.com | www.stablelab.com

¹ Oertly, et al. The accuracy of Serum Amyloid A in determining early inflammation in horses following long-distance transportation by air. AAEP Proceedings, 2017 460-461.

All trademarks are the property of Zoetis Services LLC or a related company or a licensor unless otherwise noted. Stablelab is a registered trademark of Epona Biotech Limited, used under license. © 2020 Zoetis Services LLC. All rights reserved. STB-00074

